

NF- κ B reporter assay

1. plate cells in a 6-well plate.
(100% confluent 10cm plate = 14wells of 70% confluent 6well plate if you grow it for o/n)
2. mix DNA (including GFP 20ng, renilla 5ng, firefly 10ng up to 3ug – use smaller total volume)
3. transfect cells with 3ug DNA and 7ul of lipofectamine per well
 - ① dilute 3ug of DNA into 250ul DMEM
 - ② dilute 7ul of lipofectamine into 250ul DMEM incubate for 5' @ r/t
 - ③ mix DNA and lipofectamine , incubate for 20' @ r/t
 - ④ add 500ul of mixture per well
 - ⑤ change medium after 4-6hrs (optional)
4. harvest cells
 - medium suction -> wash with PBS -> add 1ml PBS and blow off cells -> collect cells in e-tube
 - > spin @ 1200rpm for 1' -> pour off supernatant-> add 500ul PBS-> split to two tubes
 - > spin down cells-> pour of supernatant -> freeze at -20C
- 5-1. Dual luciferase assay
 - ① add 100ul of 1x PLB(passive lysis buffer)/tube
 - ② vortex it gently, leave for 10'
 - ③ spin at max speed for 1' (supernatant- cytosol, pellet – nucleus)
 - ④ use 20ul for each assay
 - ⑤ make Stop& Glo reagent (prepare it fresh)
 - ⑥ check luciferase activity
 - turn power on -> others-> oper.func->reagent-> prime *2
 - measure -> protocols-> 5-> enter->yes
 - when ended-> exit-> others-> oper.func-> reagent-> prime *2 -> prime with dH2O *2 -> prime with air*2 -> turn off
- 5-2 Quantitative Western
 - ① cell lysis in 50ul of 1xNET, 0.5% Triton X-100, 1mM PMSF at r/t
 - ② spin at max speed for 2' , take supernatant
 - ③ add sample buffer and boil at 95-98C for 2'
 - ④ run at 9% SDS-PAGE for running, 5% SDS-PAGE for stacking gel
 - ⑤ transfer to nitrocellulose memb
 - ⑥ block with 5% milk
 - ⑦ primary Ab (1:500) for 1hr
 - ⑧ secondary Ab (1:1000) for 1hr
6. Alkaline Phosphatase reaction
 - add 66ul of NBT stock + 66ul of BCIP to 10ml AP buffer
 - develop at r/t for at least 10min
 - to stop, wash the memb under tap water
7. Quantitate the amount of protein by using Image Quant software